


# Phasor™


## 3D Holography Photomanipulation Laser System




**Simultaneous photostimulation of multiple ROIs with complex shapes**



**Photostimulation of ROIs in 3D**




**Photostimulation of ROIs in 3D at different z positions**



- ✓ Complex ROIs with nm resolution
- ✓ Simultaneous photostimulation of multiple ROIs
- ✓ High power density: all incident laser power sent to SLM
- ✓ Fast update rate (~60Hz)

**3D capability:**

- ✓ Photostimulation at different planes than the focal plane
- ✓ Simultaneous photostimulation of ROIs in different z positions
- ✓ Simultaneous imaging and photostimulation at different planes



With diffraction limited spots